

The Safety Sleeper™

1-724-96-SLEEP
TheSafetySleeper.com

MODEL OPTIONS

Twin or Full sizes available

		Base 100	Standard 200	Deluxe 300	Professional 400	Premium 500
	ONE STANDARD SIDE ENTRY WITH LOCKING SYSTEM Picture shows the standard entry option of The Safety Sleeper™ with the door fully opened. *	✓	✓	✓		
	CONVERTIBLE SIDE ENTRY WITH TWO ZIPPER ACCESS AND LOCKING SYSTEM ON BOTH SIDES. This model gives the most access options for the caregiver and allows for access on both sides and above the bed. *				✓	✓
	END OF BED ENTRY: Additional entry at the foot of the bed allows for a second access point. The side entry is installed to the left when facing the end of bed entry. *			✓		✓
	ADDITIONAL COVERLET: All models of The Safety Sleeper comes with one coverlet which functions as a zip on sheet. Models with an additional coverlet ensure easy transitions between coverlet washes. ** Waterproof mattress pads are also available for purchase as the coverlets are not waterproof.		✓	✓	✓	✓
	FRAME PADS: (bar cushions) cover the top and sides of the aluminum frame and add additional safety and comfort to your bed. Recommended for seizure diagnoses or self-harming behaviors.		✓	✓	✓	✓
	STABILIZING STRAPS: Straps keep bed from tipping or moving, recommended for older children, adults or to counteract rambunctious activity. PLEASE NOTE: Stabilizing straps need to be attached to your traditional bed frame or box spring.	✓	✓	✓	✓	✓
	ACCESS POINTS: 1-1/8" ID metal grommets installed on left and right sides of the enclosure to serve as access points for use with feeding tubes, oxygen tubing, pulse-ox sensors, etc. Two on each side, four total.			✓	✓	✓
	ATTACHMENT LOOPS: The six loop attachments are installed on the roof of the enclosure to be used for items intended for hanging.					✓
	INSIDE POCKET: 8" X 10" in size and made from black mesh fabric, the pocket can be used for storing communication devices, toys, eye glasses or other items.					✓
	AIR MATTRESS FOR TRAVELING: Comes with built-in electric pump for US & Canada and battery pump for other international sizes. ***	✓	✓	✓	✓	✓
	TRAVEL CASE: Our beds weigh less than 38 lbs / 17 kg and come with a durable travel suitcase. ***	✓	✓	✓	✓	✓

* All entries come with zipper access and locking system

** Additional coverlets sold separately

*** Not covered under warranty

The Safety Sleeper
is a product of



The Safety Sleeper™

A Safe, Fully Enclosed and Versatile Bed System For Children And Adults With Special Needs

Q: Who was The Safety Sleeper developed for?

A: Anyone living with a diagnosis that affects sleep or having a night-time sleep schedule. These individuals may experience wandering at night, elopement, falling out of bed during seizures, night tantrums, head banging, inability to self-soothe due to over-stimulation; hurting self, damaging the home or causing fires due to inability to discern safety issues.

Q: What are some diagnoses that benefit from The Safety Sleeper?

A: Alzheimer's Disease, Angelman's Syndrome, Autism, Brain Injury, Cerebral Palsy, Dementia, Down Syndrome, Epilepsy, Genetic Disorder, Mitochondrial Disease, Sleepwalking, Smith Magenis Syndrome (SMS).

Q: Is The Safety Sleeper covered by insurance?

A: In many states, The Safety Sleeper has been covered through Medicaid, Medicaid HMOs, Medical Waiver and in several cases – private insurance. Each state is different in coverage and most insurers require a physician's prescription and letter of medical necessity. Our Customer Care team is available to assist you with resources to help with the insurance approval process. Call us at 724-96-SLEEP or email sales@thesafetysleeper.com.

Q: How does The Safety Sleeper compare in price and value?

A: Less than \$4,000, The Safety Sleeper is the most affordable bed on the market against comparable bed systems. Lower priced options exist, but they do not offer long lasting durability and protection and are meant for only small children. Other more comparable beds range in cost from \$4,000 – \$12,000.

Q: How easy is The Safety Sleeper to assemble?

A: The Safety Sleeper should take 15 minutes to assemble. Instructions are included with the product and a link to our assembly video is provided. Here are two important tips about the assembly of The Safety Sleeper: 1.) The Coverlet functions as the sheet for your bed. It is designed to fully enclose the mattress and has a hidden zipper. There is a separate video about the coverlet on our website. 2.) The air mattress is for TRAVEL USE ONLY. For daily use, a standard twin/full size mattress between 9-11 inches high is recommended.

FDA Registered: Registration # 3011389091 | Owner/Operator # 10048434. As a Class 1 medical device, The Safety Sleeper is designed to meet the FDA 7 Zones of Entrapment guidelines of the most common entrapment areas.



Different Color Options Available



Bed material is Fire Retardant